

Frequently Asked Questions



Are the babies comfortable on the mattress covers?

Yes the material used is similar to the Martex washable incontinence pads, which means any sweat or fluid from the baby is absorbed into the cover, rendering the baby dry.

Are the colours too dark?

No, the parents of the babies like the patterns, which look less clinical than the white sheets. The colours absorb the ultra violet light during phototherapy, which is better than reflecting back onto the Nurse.

Are the mattress covers washable?

Yes, using the same washing instructions as for the Martex Bed Pad

SCP Clamps: How do we fit them?

Place the endotracheal tube through the clamp, assessing the amount that needs intubating, cut off the tube and fit the connector. Intubate the baby and then tighten the clamp. Fix the tube and clamp in position by tying to the bonnet

How do we loosen the clamp?

Use Spencer Wells forceps and loosen the clamp by pulling the left side of the serrated part, away from the other side of the serration