

Frequently Asked Questions



What pressures can be generated by the pump?

Whatever water pressure you need to irrigate the ureter or kidney. This is dependent upon the amount of pumping done with the foot. Always ensure the outlet is free for draining

How can I purchase the Peditrol®

Either purchase the pump and purchase the disposables, alternatively take the pump on permanent loan, with the proviso that you keep purchasing the disposables

Does the product need to be sterile?

The disposables are double wrapped and sterile and the Pump has to be clinically clean.

Can this be used in conjunction with a Laser Foot Pump?

Yes and with any other kind of Foot pump. It is advisable to use the opposite foot to stop any confusion.